

Invia® Lite™

Negative Pressure Wound Therapy




“ Safe, simple, **effective**, accessible
Negative Pressure Wound Therapy for
your **Long Term Care** patients from the
brand you know and trust ... **Medela.** ”



Why Negative Pressure Wound Therapy?

- Promotes wound retraction
- Provides mechanical stress
- Removes exudate and infectious materials
- Reduces edema
- Promotes vascular perfusion and blood flow
- Promotes granulation tissue formation by facilitating cell migration and proliferation
- Provides moist wound healing environment



Invia® Lite™ Wound Therapy...

Clinical efficacy

- Negative Pressure Wound Therapy backed by clinical research*
- Non adherent wound interface allows distribution of negative pressure without tissue integration
- Kerlix™ AMD™ Antimicrobial gauze for reducing bacterial growth in pooling layer

Safe and reliable pump

- 50 years of experience in Medical Vacuum Technology
- Measures vacuum at the drain site and at the pump assuring prescribed therapy
- Alarm system
- Easy dressing changes may reduce nursing time



Medela Invia® Lite™

Safe, simple, **effective**, and accessible

Negative Pressure Wound Therapy for your
Long Term Care patients featuring quality Swiss
technology from a brand you know and **trust**

Requirements for State-of-the-Art NPWT	Product Characteristics	Benefits
Clinical Validation	<ul style="list-style-type: none">• Pump backed by 50 years of Swiss vacuum technology• Clinically validated negative pressure system	<ul style="list-style-type: none">• Proven reliability• Superior pump technology• Positive clinical outcomes
Safe and Reliable	<ul style="list-style-type: none">• Alarm system• Patient lock-out mode• Double lumen tubing quickly detects clogs or kinks	<ul style="list-style-type: none">• Early detection and intervention of problems or irregularities in the course of therapy• Digitally monitors therapy
Comfort	<ul style="list-style-type: none">• Non-adherent wound contact layer• Antimicrobial gauze• No-sting barrier film wipes• Strip paste• Dressing change 2-3 times weekly	<ul style="list-style-type: none">• May minimize trauma to new tissue• May reduce the risk of infection• May reduce patient discomfort during dressing changes
Mobility	<ul style="list-style-type: none">• Complete system with integrated suction source• Lightweight 2.2 lb. pump• Discreet shoulder carrying case• Universal rail holder	<ul style="list-style-type: none">• Greater patient mobility and freedom• Assists patient in returning to activities of daily living
Ease of Use	<ul style="list-style-type: none">• Digital screen• Push button vacuum control• Intuitive menu controls	<ul style="list-style-type: none">• Clinician can easily make changes in therapy• Patient cannot accidentally change therapy

Invia® Lite™ Wound Therapy helps promote the healing of:

Diabetic/Neuropathic Ulcers • Pressure Ulcers • Chronic Wounds • Acute Wounds • Dehisced Wounds • Grafts

