

**NPD 1000™**  
NEGATIVE PRESSURE WOUND THERAPY SYSTEM



**CONVENIENCE**  
that produces **POSITIVE OUTCOMES**

**Kalypto**  
Medical

## THE CONVENIENT WAY TO ADDRESS WOUNDS THAT CAN OTHERWISE TAKE A TOLL.

Lower extremity ulcers are the most common chronic wounds with 3.2 million Americans affected at a cost of 2 million lost work days annually and greater than \$1 billion in prescriptions.<sup>1</sup> Kalypso Medical offers a convenient negative pressure wound therapy (NPWT) system to help heal chronic wounds, while helping patients maintain their lifestyles.



### A LIGHTWEIGHT PUMP THAT'S PROGRAMMED FOR CONVENIENCE.


- 3-button programming for the clinician
- 1-button operation for the patient
- Visual and audible alarms
- 3+ weeks of power with 3 AA batteries
- Features continuous or intermittent negative pressure modes

# EASY-TO-USE ANSWER FOR

## ONE-PIECE, PEEL-AND-STICK DRESSING DESIGNED TO ACCELERATE APPLICATION TIME AND WOUND HEALING POTENTIAL.



*Applies in just minutes.*




*Easy, self-sealing gasket*

*Nonwoven absorbent pad  
"wicks and locks" exudate,  
eliminates the need for a canister*

*Silverlon® wound contact layer  
(minimally adherent  
antimicrobial barrier)*

**ADVANCED VERSATILITY  
TO EXPAND YOUR CLINICAL  
TREATMENT OPTIONS.**

- Works on a variety of wounds, including those not typically treated with conventional NPWT:
  - Diabetic foot ulcers
  - Venous insufficiency ulcers
  - Pressure ulcers
  - Arterial insufficiency ulcers
- Use with compression therapy or over grafts and flaps



*The NPD 1000 is intended to remove a small amount of fluid from wounds, which may promote healing.*



*Treat adjacent wounds simultaneously.*



*Compression therapy-compatible.*

# HARD-TO-HEAL WOUNDS.



**LIGHT AND COMPACT ENOUGH TO  
HELP ALL PATIENTS MAINTAIN THEIR  
NORMAL LEVEL OF ACTIVITY.**










The NPD 1000 system reduces the burden of negative pressure wound therapy.

- Pump is lightweight, at 8 ounces
- Pump is compact – 3.2" w x 5.1" h x 1.3" d
- Discreet due to its size and quiet nature



# ORDERING INFORMATION

## NPD 1000™ NEGATIVE PRESSURE WOUND THERAPY SYSTEM

| Product   | Order Code   | Description   | Size                               | #/Case |
|---|--------------|---|------------------------------------|--------|
|  <b>25cc Oval Dressing</b>         | M1-05X08-S05 | 25cc fluid capacity   | 2.1" x 3.2" pad                    | 5      |
|  <b>50cc Oval Dressing</b>         | B1-08X12-S10 | 50cc fluid capacity   | 3.0" x 4.6" pad                    | 10     |
|  <b>50cc Forefoot Dressing</b>     | F1-11X16-S10 | 50cc fluid capacity   | 3.0" x 4.0" pad<br>1.3" x 9.0" pad | 10     |
|  <b>140cc Dressing</b>             | B1-12X16-S10 | 140cc fluid capacity  | 4.8" x 6.3" pad                    | 10     |
|  <b>70cc Sacral Dressing</b>       | S1-11X12-S05 | 70cc fluid capacity   | 4.8" x 4.5" pad                    | 5      |
|  <b>NPD 1000 NPWT System Pump</b> | P-1000i-K    | Pump module, control module, 3 AA batteries, tubing w/Quick Connector | 3.2" w x 5.1" h x 1.3" d           | 1      |
|  <b>Y-Connector</b>              | T-4000-05    | Attaches 2 dressings to 1 pump  |                                    | 5      |
|  <b>Carrying Case</b>            | PC-1000-1    | Designed to fit the NPD 1000, with belt loop, carrying strap          |                                    | 1      |
|  <b>Dressing Tubing - Long</b>   | T-2000-05    | Attaches dressing to pump tubing                                      | 18"                                | 5      |

**Rx only.**

Contact your local Kalypto Medical Representative today for more information on the NPD 1000™ NPWT System.

**Reference:**

1. <http://www.centerforskinintegrity.com/html/educating.htm>. October 22, 2010.

