


An Advanced Sacral Dressing for Negative Pressure Wound Therapy (NPWT)



Sacral 70cc NPWT Dressing Specifications:

Dressing size: 8.7" x 8.6" (22.4cm x 22cm)

Wound Contact Area: 4.47" x 4.73" (11cm x 12cm)

Kalypto Medical NPD 1000 NPWT System

Easy. Portable. Effective.

Effectively use the Kalypto Medical NPD 1000 NPWT System for a variety of full-thickness wounds and patient types





Dressing:

- Easy-to-Use, advanced dressing technology
- Rapid application
 - Multiple dressing options for various wound types
 - Wicks and locks exudates

Pump:

- Portable, lightweight and ultra-quiet
- May increase patient compliance
 - Powered 3+ weeks by 3 AA batteries

Product	Part Number	Name	Description
	S1-11x12-S05	NPD 1000 70cc Sacral Dressing	70cc Fluid Capacity Sacral Dressing Size: 8.7"x 8.6" (22.4cm x 22cm) Wound Contact Area: 4.47"x4.73" (11 cm x 12cm) Pigtail Length: 18" CASE OF 5
	P-1000i-K	NPD 1000 System Device	NPD 1000 System Device (Pump & Control Module) Dimensions: 3.2" W x 5.1" H x 1.3" D; Weight: 8 oz. Device Tubing Length: 18" 3 AA Batteries

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician

Contact your local Kalypto Medical Representative today for addition information on the NPD 1000 NPWT System.

